

## UM-AJV5-01-1060

### Revision history

Rev.	Date	Edited by	Sections affected	Approved by
0.10	2007.06.27.	T. Józsi, L. Süle	Preliminary version	
1.00	2007.07.05.	C. Magyar	Revision	
1.03	2009.03.23.	C. Magyar	power supply models	
1.04	2009.10.13.	I. Losonczy	New pump, new power supply	
1.05	2012.02.24.	C. Magyar	cover page, p7:Front panel, p9-10 AXIOHM printer added p11: power supply image changed	
1.06	2012.12.3	M. Switzer	Update Address, 5-part and 3-part Species	

	Date	Name	Signature
Issued	12/3/2012	Mike Switzer	
Checked and approved		Árpád Gyetvai	

# *Abacus junior vet5*

## Hematology Analyzer


### User's Manual




1.06 Release



DIATRON MI PLC  
H-1097 Budapest,  
Tablas utca 39. HUNGARY  
Tel.: +36 1 43 9800  
Fax: +36 1 436 9809  
[www.diatron.com](http://www.diatron.com)  
[support@diatron.com](mailto:support@diatron.com)

## Contents

<b>1. INTRODUCTION .....</b>	<b>4</b>
1.1 Intended Use .....	4
1.2 The Instrument .....	4
1.2.1 Patient Testing .....	5
1.2.2 Reagents .....	6
1.2.3 Technical Operation .....	6
1.2.4 Calibration and Quality Control .....	6
1.3 Instrument features .....	7
1.4 Parts of the Analyzer .....	11
1.4.1 Function of the Fluidics .....	11
1.4.2 Control Panels .....	12
1.4.3 Display .....	12
1.4.4 Keyboard .....	13
1.5 Control Material .....	15
1.6 Accessories .....	15
1.7 VET5 Specifications .....	16
<b>2. INSTALLATION .....</b>	<b>18</b>
2.1 General information .....	18
2.2 Environmental factors .....	18
2.2.1 Electrical requirements .....	18
2.2.2 Space requirements .....	19
2.2.3 Peripherals .....	19
2.2.4 Reagents and waste handling .....	20
2.2.5 Maintenance .....	20
2.2.6 Cleaning .....	20
2.2.7 General points .....	20
2.3 Unpacking and installation .....	21
2.3.1 Turning the Instrument ON .....	23
2.3.2 Turning the Instrument OFF  .....	24
2.3.3 Preparing for shipment .....	25
2.3.4 Emergency handling .....	25
<b>3. MENU SYSTEM .....</b>	<b>26</b>
3.1 General Information .....	26
3.1.1 Navigating in the Menu System .....	26
3.1.2 Menu structure .....	28
<b>4. OPERATING PRINCIPLES .....</b>	<b>33</b>

4.1	Impedance Method .....	33
4.2	Principle of HGB Measurement .....	33
4.3	Parameters.....	34
4.4	Absolute and Linearity Ranges of Measured Parameters .....	35
<b>5.</b>	<b>ROUTINE UTILIZATION and MEASURE</b>  .....	<b>36</b>
5.1	Sample handling .....	36
5.2	Sample analysis.....	39
5.2.1	Sample preparation .....	39
5.2.2	Modifying lyse quantity .....	39
5.2.3	Sample information .....	40
5.2.4	Results .....	40
5.3	Measure.....	43
5.3.1	Measure local menu .....	43
5.3.2	Blank Measurement.....	43
5.3.3	Using Pre-diluted Mode .....	45
<b>6.</b>	<b>DATABASE</b>  .....	<b>47</b>
<b>7.</b>	<b>UTILITIES</b>  .....	<b>49</b>
7.1	Maintenance.....	49
7.1.1	Regular Maintenance Jobs.....	49
7.1.2	Weekly maintenance.....	51
7.2	Calibration .....	55
7.3	Quality Control Procedure .....	57
7.4	Diagnostics.....	59
7.5	Settings.....	59
7.5.1	Printer settings.....	60
7.5.2	Customize.....	64
7.5.3	Units .....	66
7.5.4	Date and time .....	68
<b>8.</b>	<b>Printing</b> .....	<b>69</b>
<b>9.</b>	<b>Fluidic Schematics</b> .....	<b>72</b>
<b>10.</b>	<b>Warning labels on the analyzer</b> .....	<b>73</b>

# 1. INTRODUCTION

## 1.1 Intended Use

**Abacus junior vet 5** hematology analyzer is a fully automated cell counter designed for *in vitro* diagnostic use.

The compact instruments were developed for small clinics, point-of-cares and vet offices.

## 1.2 The Instrument

**Abacus junior vet 5** is a fully automated, bench top hematology cell counter.

It implements the so-called Coulter-method for counting cells passing through a small aperture, and measures the hemoglobin content of red blood cells.

The analyzer features a graphical LCD display module and a foil keypad of 29 keys including 6 software buttons (with icons), 6 function keys (above LCD) and has a START button.

The instrument allows sending results to an external printer, or can have an optional built-in printer module.

Its internal memory is capable of storing 2000 records with full histograms, and individual patient data. QC measurements can also be performed, and stored. The software operating the instrument is easy to upgrade using a USB drive. The instrument allows connecting to a host computer for uploading records stored in the memory through a USB SLAVE port or serial link, and also enables archiving and restoring records to and from USB drive.

---

<b>NOTE:</b>	<b>If the equipment is used in a manner different from which the manufacturer specified, the protection provided by the equipment may be impaired. Misuse of equipment or use other than its intended purpose will invalidate conditions of warranty. The accuracy and precision may also be impaired.</b>
--------------	--

---

### 1.2.1 Patient Testing

The analyzer can process 24-30 samples per hour in 3-part mode, 16-20 samples per hour in 5-part mode. Samples can have individual sample data, and additional parameters.

Results can be printed on the optional internal or external printer. The print-out format can be customized by the user.

**Abacus junior vet 5** determines 22 hematology parameters including five-part WBC differential in 5-part mode. The instrument requires 50 µl (2 x 25µl) of the whole blood sample:

WBC	total white blood cell count
LYM	lymphocytes count
MON	monocytes count
NEU	neutrophil count
EOS	eosinophil count
BAS	basophil count
LYM%	lymphocyte percentage
MON%	monocytes percentage
NEU%	neutrophil percentage
EOS%	eosinophil percentage
BAS%	basophil percentage
HGB	hemoglobin
RBC	red blood cell count
HCT	hematocrit
MCV	mean corpuscular volume
MCH	mean corpuscular hemoglobin
MCHC	mean corpuscular hemoglobin concentration
RDWc	red cell distribution width
PLT	platelet count
PCT	platelet percentage
MPV	mean platelet volume
PDWc	platelet distribution width

### 1.2.2 Reagents

Only reagents supplied by the manufacturer should be used with the analyzer, otherwise accuracy cannot be guaranteed.

<b>Diluent:</b>	Isotonic saline solution used to dilute whole blood specimens and to rinse the fluidic system between measuring procedures.
<b>Lysing reagent:</b>	Used to create hemolysate for 3-part WBC differential and for total WBC and HGB.
<b>Cleaner:</b>	Used to perform cleaning process of the fluidics.
<b>Rinse:</b>	Used in cleaning process of the fluidics.
<b>VET5:</b>	Used to dilute whole blood and differentially hemolyse white blood cells to separate eosinophil granulocytes from other WBC by volume. It is suitable to determine EOS, EOS%, BAS and BAS% parameters.

### 1.2.3 Technical Operation

As **Abacus junior vet 5** is a fully automated instrument, operating it requires minimal training or technical support. Operator interaction is reduced to the following:

- Perform a *Blank Measurement* in case the instrument is not used for a specific time
- Enter sample and/or patient data
- Apply the sample for analysis
- Print results either one-by-one, or in groups by selecting records from the database
- Perform simple weekly maintenance, as described later in this description (7.1).


### 1.2.4 Calibration and Quality Control

**Abacus junior vet 5** arrives at your laboratory factory-calibrated and ready to use. However, calibration needs updating whenever you find that the results have slightly changed, or a different or new control material is used. With each control material you receive to be used with the instrument, you will find a control sheet listing the parameters the instrument should match. Perform these calibrations as explained in a later chapter (7.2).

Quality Control is used for checking for proper calibration and performance of the analyzer. Running these samples should happen on a regular basis, as also explained in a later chapter (7.3).

## 1.3 Instrument features

Figures 1 and 2 show front and rear view of **Abacus junior vet 5** hematology analyzer.



**Figure 1. Front view**

1. USB connector
2. OK key
3. Numerical keypad
4. Function keys
5. Graphic liquid crystal display
6. HELP key
7. Measure function key
8. Database function key
9. Utilities menu key
10. Printing function key
11. Exit menu key
12. Cursor control keys
13. Status indicator
14. START key
15. Sample rotor (with interchangeable adapters)





**Figure 2. Rear view**


1. Reagent tubing connections
2. On/Off switch
3. External power supply inlet 12VDC
4. USB A port
5. Serial port (RS 232)
6. USB B port
7. Keyboard connector (PS2)

The analyzer is equipped with built-in printer allowing immediate report generation. Abacus Juinor VET5 is designed to support two printer models; a Seiko and an Axiohm model. Both printers use 58mm wide thermal roll paper.

Installation of the printer is a factory option. Each analyzer that has the built-in printer option will be equipped either with the SEIKO or with the AXIOHM printer module.

The SEIKO model has a flat top panel with color matching that of the instrument cover. To open its lid, press the black button. Just let the thermal printer paper fall into the housing with its open end toward the black roller. Holding the open end of the paper, close the lid so that the end of the paper gets clipped between the black paper guide and the printer mechanics.

You can select to use this printer in the “**Printer Settings**” menu (see chapter 7.5.1).



**Figure 3. Built-in printer**

The AXIOHM printer has a black cover. It also has a blue status light, that blinks when the printer paper is out.

**To replace paper in the printer:**

- open the paper lid (pull the lid upwards by the handle)
- remove central plastic roller of old paper roll
- unwind new paper roll, so that the “starting edge” is coming from down under towards you



- gently drop the new roll into the holder of the printer, and hold the “starting edge” with your hand, and make sure it comes out on the front of the printer



- close the lid, making sure that the paper is captured between the lid and the front of the printer



The analyzer works with an external power supply. The next figure shows the power supply unit generating 12VDC. The power supply module has a so-called auto range input, allowing operation on 230V or 115V power system. The power supply complies CE and UL safety certifications.



The input socket is a standard power cable connection, the output is a special, lockable socket as shown in the picture.

Only the provided power supply shall be used with the instrument. ("Protek Electronics Corp." Model ID: PUP80-12-N8, or "GLOBTEK INC", Model ID: GT-81081-6012-T3")

## 1.4 Parts of the Analyzer

**Abacus junior vet 5** hematology analyzer is composed of three main units:

**Fluidic System:** Performs sampling, diluting, mixing, and lysing and Lyse2 lysing functions. Generates the regulated vacuum used for moving cells through the aperture during the counting process.

**Data Processing System:** Counts, measures and calculates blood parameters, generates and stores numerical results and histograms.

**Control Panel:** Features an LCD display, a 29-button keypad, and USB and serial (computer) interfaces.

### 1.4.1 Function of the Fluidics

For the Schematics of the fluidics system, see Section 9.

#### Sample aspiration and dilution:

Stages of the blood testing process	
a.	25 µl of anti-coagulated (EDTA) whole blood sample is aspirated through the sampling needle, mixed with VET5 reagent and diluent for the 5-part differentiation.
b.	After the first counting and washing process, another 25 µl of anti-coagulated (EDTA) whole blood sample is aspirated again and mixed with 4 ml of <b>diluent</b> and stored in the chamber (mix dilution).
c.	25 µl of the mix dilution is aspirated into and stored in the needle during WBC

	measurement and the hemoglobin analysis
d.	<b>Lysing reagent</b> is added to the mix dilution held in the chamber for WBC differential analysis. This amount of <b>lysing reagent</b> is species dependent and may be changed by the operator within “Limits” menu.
e.	After WBC counting and a washing process, 5 ml of diluent is added to the second dilution (using the 25 µl of mix dilution stored in the needle).
f.	This portion is analyzed for RBC count, PLT count and their parameters.

**Table 1.**

**Dilution rates used within *Abacus junior vet 5*:**

- DIFF dilution 1:160
- Mix dilution 1:160
- RBC dilution 1:32 000
- WBC dilution 1:196 (depends on lyse amount)

**Measurement times:**

- DIFF measurement 8 seconds
- WBC count 8 seconds
- HGB measurement 3 seconds
- RBC/PLT count 8 seconds

## 1.4.2 Control Panels

### START button

Pressing and releasing the START button triggers an analysis cycle.

### Status indicator

A three-color LED is located near the START button.

Its actual color indicates the status of the analyzer.

LED color	Analyzer status
● Green	The analyzer is ready to work. Analysis can be initiated.
● Red blinking	Blood sample can be removed when the LED blinks red 3 times and the instrument beeps 3 times.
● Red	The analyzer is currently performing an analysis. No new measurement can be started.
● Yellow	The analyzer is performing a maintenance process.
● Yellow blinking	The instrument is in stand-by and display light is off.

## 1.4.3 Display

The display is a 240 x 128 dots, high contrast backlit graphic LCD (Liquid Crystal Display) module.
















#### 1.4.4 Keyboard

The foil keypad is composed of the following (shown in Section 1.3 - Figure 1):







- **Numeric keys** for entering numerical data, and selecting menu items
- **Function keys** with specific functions. These functions are menu-dependent and are indicated by icons appearing above the keys
- **Hardware function keys (short-cut keys)** for easier navigation between menus
- **Cursor control keys** ↑ and ↓ for moving between database items, ← and →, for moving between parameter columns or menu levels
- **START key** for initiating an analysis cycle
- **OK key** for confirming data
- **Del key** for deleting characters
- **Help key** for HELP function

## Function keys

Below is a list of possible icons and functions assigned to the so-called soft-keys (function keys)

Function key	Action triggered
	Exit from actual menu or action
	Leave data-entry menu without saving any changes made to it (Cancel)
	Confirm the results or changes made (OK)
	Redo action (e.g. Blank measurement)
	Display histograms of the highlighted patient ID or QC Lot No.
	Show data in table format
	Enter/modify sample/patient data
	Select between result pages
	PAGE-UP / PAGE-DOWN key in a multi-page menu
	Change scaling of Levey-Jennings chart (16 or 64 days)
	Patient type selection
	Confirm error
	Go to local menu (database, measurement)
	Limits
	Stop a running process

Below we list the icons and functions assigned to so-called hardware function buttons

Function key	Action triggered
	Information
	Measuring process at once
	Database
	Utilities menu
	Printing function
	Exit menu

## 1.5 Control Material

**Abacus junior vet 5** allows continuous monitoring with DiatroCont5 control (control blood). This should be matched to the types of samples usually run on the instrument. Specification for this material (assay values and allowed tolerances along with expiry date) is always enclosed with the approved control material.

## 1.6 Accessories

Below is a list of accessories shipped with your **Abacus junior vet 5** instrument. This list can also be referred as the “**Abacus junior vet 5 - pack**”

**Abacus junior vet 5** Hematology Analyzer

**Abacus junior vet 5** User's Manual (this booklet)

**Abacus junior vet 5** Reagent Tubes (with colored connector caps)

Diluent tube (green)

Lyse tube (yellow)

L\_VET5 tube (orange yellow)

Clean tube (blue)

Waste tube (red)

Rinse tube (white)

**Abacus junior vet 5** Cleaning Tube Kit.

**Abacus junior vet 5** Caps for reagent containers (matching reagent tube connector cap colors).

**Abacus junior vet 5** Waste Container (20 L).

**Abacus junior vet 5** External power supply and Power cable.

**Abacus junior vet 5** Tube adapters.

**Abacus junior vet 5** Thermal roll-paper.

Reagent Tubing Kit:



Cleaning Tube Kit:





## 1.7 VET5 Specifications

<b>Chambers</b>	25 µl of whole blood in 3-part mode, 50 µl of whole blood in 5-part mode
<b>Chambers</b>	50 µl of whole blood in pre-diluted 3-part mode
<b>Reagent system</b>	1 unified chamber for diluting whole blood and counting
<b>Aperture diameter</b>	Isotonic Diluent, Lyse, VET5 Lyse, Cleaner and Rinse
<b>Throughput</b>	80 µm
	22-28 tests/hour in 3-part mode, 15-20 tests/hour in 5-part mode

Characteristics	Accuracy	Reproducibility (CV)	Correlation coefficient(R)	Carry-over sample to sample	Test range	Unit
<b>WBC</b>	5%	<10%	>0.85	<1%	4.0-20.0	10 <sup>3</sup> /µl
<b>RBC</b>	5%	<10%	>0.85	<1%	4.0-15.0	10 <sup>6</sup> /µl
<b>HCT</b>	5%	<10%	>0.85	<1%	25.0-50.0	%
<b>HGB</b>	5%	<10%	>0.85	<1%	6-16	g/dl
<b>PLT</b>	10%	<20%	>0.75	<3%	200-900	10 <sup>3</sup> /µl
<b>LYM%</b>	20%	<20%	>0.75	N/A	15-50	%
<b>NEU%</b>	20%	<20%	>0.75	N/A	25-80	%
<b>EOS</b>	20%	<20%	>0.75	<0.1	0.5-2.0	10 <sup>3</sup> /µl

<b>Sampling method</b>	Open tube system with automatic sample rotor.
<b>Sample types</b>	Dog, Horse, Cow, Cat Calibrator/QC material (5-part mode) Ferret, Goat, Guinea Pig, Mouse, Pig, Primate, Rabbit, Rat, Sheep (3-part mode)
<b>Clog prevention</b>	High-voltage pulse on aperture in each analysis cycle, chemical cleaning and high pressure back-flush of the aperture using Cleaner reagent.
<b>Cleaning procedure</b>	High-voltage burst of the aperture, high-pressure back-flush, chemical cleaning of the aperture.
<b>Quality control</b>	6 QC levels, QC parameters include: mean, ± range, SD and CV for all measured and calculated parameters, 16- and 64-day Levey-Jennings charts, separate QC database.
<b>Calibration</b>	3-measurement automatic and manual calibration of WBC, HGB, RBC, PLT, MCV, RDW, MPV and EOS absolute. Independent calibration of pre-diluted mode.
<b>Multi-user feature (advanced feature)</b>	3-level multi-user operation with selective privilege levels, user identification with ID and password.
<b>User interface</b>	Easy-to-use, menu driven user interface with 6 software buttons (with icons), 6 hardware function buttons, cursor and numeric keys.
<b>Languages available</b>	English, French, Italian, German, Spanish, Russian, Hungarian.
<b>Data capacity</b>	2000 results, with RBC, PLT, WBC 3-part and EOS histogram (5-part mode).
<b>Host computer interface</b>	USB B port or RS-232.
<b>Data back-up method</b>	USB mass storage device (PenDrive™)
<b>Software upgrade method</b>	Free DiatronLab data management software.
<b>Printer drivers</b>	via USB A port using USB mass storage device (PenDrive™).
<b>Built-in printer</b>	Canon BJC, Esc/P, Esc/P2, Esc/P raster, PCL4.
<b>Display</b>	Seiko „Easy Paper Operation” thermal printer, 56 mm wide roll paper.
<b>Keypad</b>	240x128-dots, high-contrast, backlit, monochrome (white on blue) graphics LCD (liquid crystal diode).
<b>External keyboard</b>	29 foil keys + separate START button.
<b>Power requirement</b>	Standard PS/2 compatible keyboard.
<b>Power supply unit</b>	12VDC, 6A, 72W max. operating power.
<b>Operating temperature</b>	External, auto-ranging power unit for 100-120 or 200-240 VAC, 50–60Hz.
<b>Dimensions (W x D x H)</b>	59–86 °F (15–30 °C). Optimal temperature is 77 °F (25 °C). 12.6 x 10.2 x 14.4 in (320 x 260 x 365 mm).

**Net weight**

12 kg.

## 2. INSTALLATION

### 2.1 General information

This chapter provides instructions for the installation of the **Abacus junior vet 5** hematology analyzer. The procedures described below must be followed correctly to ensure proper operation and service. Please carefully read and follow all instructions in this *User's Manual* before attempting to operate **Abacus junior vet 5**.

**Abacus junior vet 5** hematology analyzer is a precision instrument, and must be handled accordingly. Dropping or other improper handling of the instrument will disturb calibrated mechanic and electronic components and/or cause other damage. ***Always handle the instrument with care.***

### 2.2 Environmental factors

Operate **Abacus junior vet 5** within the ambient temperature range of 15-30°C and relative humidity of 65% ± 20%. The optimum operating temperature is 25°C.

Avoid using the instrument in areas of extreme high or low temperatures or where it is exposed to direct sunlight. If it is kept at a temperature less than 10°C, the instrument should be allowed to sit for an hour at the correct room temperature before use.

Reagents should be stored at a temperature range of 15-30°C.

The instrument should be placed in a well-ventilated location. The instrument should not be placed near potentially interfering devices capable of emitting radio frequencies (e.g. radio or television receiver, radars, centrifuge, X-ray devices, fans, etc.).

Operation at an altitude over 3000 meters (9000 ft) is not recommended.

Instrument is designed to be safe for transient voltages to INSTALLATION CATEGORY II and POLLUTION DEGREE 2.

Environmental and electrical requirements have been provided to insure the accuracy and precision of the instrument and maintain a high level of operational safety for lab personnel.

#### 2.2.1 Electrical requirements

**Abacus junior vet 5** comes with a power cord appropriate for your power system. Proper use of the appropriate power cord assures adequate grounding of the system.

**WARNING:** Failure to proper ground the **Abacus junior vet 5** bypasses important safety features and may result in electrical hazard.

## 2.2.2 Space requirements

It is important to install the instrument in a suitable location. A poor location can adversely affect its performance. Consider the following space requirements:

- Select a location near a power source and close to a suitable drain.
- Place the unit on a clean and level surface.
- Leave at least 0.5 m space on both sides and above the instrument to access pneumatics and (optional) built in printer. A minimum of 0.2 m must be maintained between the rear panel and the wall to allow for heat dissipation and tube clearance.
- Install the reagents in a suitable place that will make your work easy. The best place is on the ground, below the supporting desk the instrument is placed on. The pneumatics system is capable of aspirating reagents from containers being 1.0 m below the reagent inputs. Make sure the reagent tubes are not bent, broken, twisted or blocked in between the desk the instrument is on and the wall behind. Such circumstances can result in instrument operation failure.
- DO NOT PLACE the reagents above the instrument, as there can be a risk of falling and spilling.

**WARNING:** Install the unit on a table or workbench. If the unit is installed without a supporting desktop under the unit, there is a possibility that **Abacus junior vet 5** could accidentally fall.

## 2.2.3 Peripherals

External peripherals should only be connected when both the instrument and the peripheral device are powered off. Possible peripherals are:

- external printer
  - o the printer must be recommended by authorized technician
  - o the printer must be approved and listed
  - o the printer must have a CE mark
- external keyboard
  - o the external keyboard must be approved
  - o the external keyboard must have a PS2 connector, or suitable adapter
- link to host computer via serial or USB port
  - o serial (USB) link cable must be approved by technician

## 2.2.4 Reagents and waste handling

Reagents should be handled according to national or international regulations.

---

**WARNING!** Reagents may cause corrosion and skin irritation. If any of the liquids leaked onto the cover of analyzer or the furniture, it has to be wiped off immediately. In case of skin contact, the liquid has to be rinsed off with plenty of water.



---

Waste generated by the unit is biohazard material. Handling and disposal must happen according to regulations regarding reagent systems.

---

**WARNING!** Waste contains poisonous substances (because of chemical content) and human origin substances causing biohazard. These substances are representing potential danger to environment. For this reason, safe handling of the waste liquid is very important.



---

## 2.2.5 Maintenance

The user should check the following components weekly:

- bottom of washing head for salt build up – should be wiped off with a damp cloth
- tubing system – by opening the side and back doors and look for any liquid leakage. If any leakage is observed, authorized technician should be contacted.

### Caution:

The following parts must NOT be opened or serviced by the user:

- electrical supply
- electronic boards

## 2.2.6 Cleaning

The instrument and its power supply should be cleaned on the outside only, using a damp cloth with a soft detergent. DO NOT let liquids get inside the units.

## 2.2.7 General points

The manufacturer guarantees work safety reliability and general characteristics under the following conditions only:

- services and repairs are performed by an authorized technician
- the electrical system of the laboratory follows national and/or international regulations
- the system is operated according to instructions contained herein

## 2.3 Unpacking and installation

1. Carefully remove the **Abacus junior vet 5** hematology analyzer from the shipping carton. Inspect the instrument for any visible signs of damage incurred during shipping. Would you find any damage, file a claim with the carrier or your distributor immediately. Check the accessories received against the packing list. Contact Service if anything is missing.
2. **CAUTION!** Prior to initial operation, allow the instrument to reach room temperature (approx. 2 hours). Rapid temperature changes in an operational unit can lead to water condensation, which may damage electronic parts.
3. Place the instrument on a firm work surface in the designated work area, near an appropriate AC electrical outlet. The connection **MUST** be grounded.

---

**NOTE**      **Before making connections:** Make sure that all power is in “OFF” setting before connections (printer, external keyboard) are made. Carefully read all literature accompanying the instrument and its accessories. Pay particular attention to the operating procedures for the external printer.

---

4. **Keyboard and external printer**  
Attach the keyboard cable to the round “KEYBOARD” port on the back of the instrument. Attach both ends of the printer cable to the appropriate ports on the printer and **Abacus junior vet 5**. Attach the AC adapter to the printer (if required) and plug it into an AC outlet.
5. **Host Computer**  
The instrument has a built-in serial port that allows connection to a host computer. Results, including histograms, may be exported. Serial (USB) I/O settings can be found in Settings.  
For installation instructions, please contact Service.
6. **Power supply**  
Connect the power supply to the instrument. Attach power cord outlet to the external power supply of **Abacus junior vet 5** and plug the other end into a properly grounded AC outlet.


**Do not switch on the instrument before connecting the external power supply to the instrument and to the AC outlet, as well as before connecting an external printer or a keyboard to the instrument.**

## 7. Reagent Containers

Place the reagent containers near the instrument, to an accessible location. Do not place the containers to a higher position than that of **Abacus junior vet 5**, because would a tube come off its connector, the fluids spoil out. Use the supplied connecting tubes and special bottle caps. Be sure that the color on each tube, cap and connector in the back of the instrument match. You can for example, place the reagent containers below the desk **Abacus junior vet 5** is installed on, as the instrument has sufficient power to draw the liquids from a lower location.

All containers should be left open (do not block the small air vent hole on the special container caps) in order to provide free airflow.

(For connections, see Figure 5.)



**Figure 5. Reagent connections**

---


**WARNING!** Reagents may cause corrosion and skin irritation. If any of liquids leaked to cover of analyzer or the furniture, it has to be wiped down immediately. In case of skin contact the liquid has to be rinsed by plenty of water.

---




### 2.3.1 Turning the Instrument ON

- In case you use an external printer (for information, read manual shipped with the printer) connect it and turn it on.
- Turn the **Abacus junior vet 5** on using the power switch on the rear panel. The On position is marked by the I symbol.



During start-up, the following screen is displayed. The software version number appears few seconds later, when the software starts.

PatID	SmpID	Date	
□8313	19	12.06.2007	13:40
□5612	20	12.06.2007	14:06
□3013	21	12.06.2007	14:16
□5630	22	12.06.2007	14:21
□5627	23	12.06.2007	14:32
□5457	24	12.06.2007	14:38
 (#9999)			

An important feature of the instrument is that when the software version is completed, the DATABASE will be displayed without any pneumatic initialization (default setting). Pneumatic movement will be initiated only when necessary for the relating process.

The default setting can be changed at Service Menu level, in this case the instrument will start with pneumatic initialization giving the possibility to perform a measuring process immediately. Please call the Service Personnel in case you want to change the setting.

---

**CAUTION! Wait 5 minutes before initiating any measuring process to allow the instrument to reach the optimal working temperature.**

---

In some cases, a priming cycle is necessary prior to sample introduction. The instrument will perform the cycle automatically if the fluid sensors are on and additional liquid in the tubing system is required. A priming cycle should be run:


- at installation
- in case of extended time out of use
- after replacement of any component related to the Fluidic System
- reagents are replaced with the instrument turned on



### 2.3.2 Turning the Instrument OFF

The instrument should never be switched off by simply flipping the power button on the rear panel. Doing so may result in erroneous operation during later use. It can be so, because the instrument uses different kinds of solutions, one of which is the so-called diluent. This liquid is an isotonic saline solution containing salt. Would it not be washed out of special units of the instrument, or would chambers not be filled with this solution may lead to dust condensation or salt build up. Therefore always follow the instructions below when switching the instrument off.

Press EXIT key  on the front panel. The following screen is displayed.



From the Exit menu, select **Shut down**. When prompted, press  to confirm.

**Abacus junior vet 5** will display a message and give a tone indicating that it is safe to shut it off.

Turn **Abacus junior vet 5** off using the power switch on the rear panel. The Off position is marked by **O** symbol.

### 2.3.3 Preparing for shipment

The second item in the shut down menu should be used when the instrument is to be shipped or left unused for a longer time (more than 2 weeks). The instrument will ask you to use the cleaning tube kit and 100ml of distilled water.

Follow the instructions appearing on the display.

! Message 5001/19300

Remove reagent tubing at rear reagent inputs (Diluent, Lyse, Lyse Vet5, Rinse and Cleaner). Leave waste connected.



EXIT

#### PREPARING FOR SHIPMENT (2)

Here, the user is instructed to remove tubing connectors, so they could be drained.

Leave the waste connector attached.

! Message 5002/19300

Connect min 100 ml distilled water to reagent inputs using cleaning tube kit.



Next, you should connect the cleaning tube kit to the reagent inputs, submerging the free end in a bottle containing at least 100 ml of distilled water.



#### Cleaning tube kit connected

The analyzer will flush any remaining reagents from the system into the waste container.

! Message 5003/19300

Remove cleaning tube kit. Keep reagent inputs free. Leave waste connected



As a next step, the analyzer asks you to remove the cleaning tube kit.

When finished, the analyzer prompts you to power off the system. Remove the waste connector after shutting down.

### 2.3.4 Emergency handling

In case of emergency situation - like instrument catching on fire (short-circuit, etc.) - cut off power immediately and use a fire-extinguisher if necessary.

## 3. MENU SYSTEM

### 3.1 General Information

This chapter contains information about the structure and usage of the software implemented menu structure.

This integrated software controls instrument operations including calculation and evaluation of measured data, displaying results and information screens storage and recalling of data.

#### 3.1.1 Navigating in the Menu System


The instrument uses a menu system to initiate actions and allow access to settings.

There are four possible ways to navigate between menus and menu items:




The most important issue is that there are 6 function keys above the LCD screen. You use these buttons to navigate between main functions of the instrument.



These keys are short cut keys: by pressing any of them you can directly access main functions, whichever submenu you are in. If you use an external keyboard you can access these functions with keys F8 through F12. (described in Menu Structure as well).



By pressing the short-cut keys, you can activate main functions (Information, Measure, Database, Utilities, Printing, Exit) of the instrument. Below the short-cut keys (in the upper part of the screen) you can find tabs showing which main function is active or which are accessible directly.

- a. You can select a menu item (move highlight) with the  $\uparrow$  ☐ and  $\downarrow$  keys and press the **OK** key to enter or activate the highlighted item. Within a submenu, you can press the  function key to return to the previous menu level.  
This method is suggested while learning instrument operation.
- b. Pressing the numeric key corresponding to the desired menu item allows selection and access of an item without the need to additionally press the **OK** key.  
Pressing the **0** (zero) key has the same effect as the  function key.  
This method is more efficient when the user knows the menu structure.
- c. You can also move between the different menu levels using the  $\leftarrow$  and  $\rightarrow$  keys. These have the same effect as **OK** and , respectively.

If a menu item would open up a submenu, then the item is indicated with a  $\blacktriangleright$  symbol on the right of the menu line.

Some results can be displayed in table format. The following keys may be used for browsing the database:

<b>3</b> $\uparrow$	page up
<b>9</b> $\downarrow$	page down
<b>1</b> $\uparrow$	jump to top of list
<b>7</b> $\downarrow$	jump to bottom of list

Several menus have items with boxes in front of the text. These indicate two-state options. The *selected (active)* state is indicated with a filled box, the *deselected (inactive)* state is indicated with an empty box. Pressing **OK** toggles its state.


Other items have circles in front of the text. These are called “radio-buttons”.

They are divided into groups separated by horizontal lines.

The function of these groups is that only one item of the group can be selected, and this is indicated with a filled circle in front of the selected item. Selecting an item of the group will move the filled circle in front of this item, emptying the circle of the old selected item.

### 3.1.2 Menu structure

The following outlines the Abacus junior vet 5's menu functions.





## Database



1 – Go to specified record

Selects one record by time/date, Sample ID, and/or Patient ID (more criteria narrows search).

2 – Selection



1 – Select by date, time and ID

Selects one or more records by range of times/dates, sample IDs, Selects all records.

2 – Select all

Deselects all records.

3 – Deselect all

3 – Change sort order



1 – Unsorted

List records in the order entered (most recent to oldest).

2 – Sort by time

Lists records by measurement time (oldest to most recent).

3 – Sort by sample ID

Lists records by sample ID (highest to lowest).

4 – Sort by patient ID

Lists records by patient ID (highest to lowest).

4 – Manage selected records



1 – Send selected records

Sends selected records to a computer through a serial connection.

2 – Delete selected records

Deletes selected records from the VET5 database.

3 – Backup selected records


Saves selected records to USB drive.

5 – View external

Reads stored records from USB drive


6 – Backup one day

Backs up all results from specified date to USB drive (insert USB drive



Toggle between graph and table display of selected results.

Opens dialog for viewing/modifying patient information (sample ID and patient type cannot be changed).



Sample ID

Name

Patient ID

Birth


Species

Sex

Doctor




## Utilities





## Utilities







## Printer

Prints selected test results.

## Exit



- 1 – Shut down
- 2 – Preparing for shipment


Shuts off the VET5 (use if the unit will not be used for over 72 hours).

Drains and shuts off the VET5 (use if the unit will be unused for more than 2 weeks, or if it will be shipped).

## 4. OPERATING PRINCIPLES

### 4.1 Impedance Method

The impedance method (a.k.a. Coulter method) counts and sizes cells by detecting and measuring changes in electrical impedance when a particle in a conductive liquid passes through a small aperture.



**Figure 7. Impedance method**

Each cell passing through the aperture – there is a constant DC current flowing between the external and internal electrodes – causes some change in the impedance of the conductive blood cell suspension.

These changes are recorded as increases in the voltage between the electrodes.

The number of pulses is proportional to the number of particles. The intensity of each pulse is proportional to the volume of that particle. The volume distribution of the cells are displayed on diagrams: WBC (EOS), RBC, and PLT histograms.

### 4.2 Principle of HGB Measurement

The lysed 1:196 sample dilution can be measured by a cyanmethemoglobin method. The reagent lyses the red blood cells, which release hemoglobin.

Hemoglobin iron is converted from the ferrous ( $\text{Fe}^{2+}$ ) to the ferric ( $\text{Fe}^{3+}$ ) state to form methemoglobin, which combines with potassium cyanide (KCN) to produce the stable cyanmethemoglobin, or hemoglobincyanide. Subsequently, the HGB concentration is measured photometrically.




---

**Note:** The above-mentioned measuring method is used to determine the HGB concentration. The HGB concentration can be measured using cyan-free lysing reagents as well. In this case the effect is the same but the used lyse is environmental-friendly reagent.

---

## 4.3 Parameters

**Abacus junior vet 5** measures and calculates 22 parameters, listed below. For each parameter we list the name, abbreviation and measurement unit in the first column. Short description for each parameter is in the second column.

White Blood Cells – <b>WBC</b> (cells/l, cells/ $\mu$ l)	Number of leukocytes $WBC = WBC_{cal} \times (\text{cells/l, cells}/\mu\text{l})$
Red Blood Cells – <b>RBC</b> (cells/l, cells/ $\mu$ l)	Number of erythrocytes $RBC = RBC_{cal} \times (\text{cells/l, cells}/\mu\text{l})$
Hemoglobin concentration - <b>HGB</b> (g/dl, g/l, mmol/l)	Measured photometrically at 540 nm; in each cycle blank measurement is performed on diluent $HGB = HGB_{cal} \times (HGB_{measured} - HGB_{blank})$
Mean Corpuscular Volume - <b>MCV</b> (fl)	Average volume of individual erythrocytes derived from the RBC histogram.
Hematocrit – <b>HCT</b> (percentage, absolute)	Calculated from the RBC and MCV values. $HCT_{percentage} = RBC \times MCV \times 100$ $HCT_{absolute} = RBC \times MCV$
Mean Corpuscular Hemoglobin – <b>MCH</b> (pg, fmol)	Average hemoglobin content of erythrocytes, calculated from RBC and HGB values. $MCH = HGB / RBC$
Mean Corpuscular Hemoglobin Concentration – <b>MCHC</b> (g/dl, g/l, mmol/l)	Calculated from the HGB and HCT values. $MCHC = HGB / HCT_{absolute}$ Unit of measurement is displayed according to the one chosen for HGB result (g/dl, g/l or mmol/l)
Red Cell Distribution Width – <b>RDW-SD</b> (fl) Platelet Distribution Width – <b>PDW-SD</b> (fl)  Red cell Distribution Width – <b>RDW-CV</b> (absolute) Platelet Distribution Width – <b>PDW-CV</b> (absolute)	The distribution width of the erythrocyte or platelet population derived from the histogram at 20% of peak  $x\text{DW-SD} = RDW_{cal} \times (P2 - P1) \text{ (fl)}$ $x\text{DW-CV} = RDW_{cal} \times 0.56 \times (P2 - P1) / (P2 + P1)$ by the factor of 0.56 CV is corrected to the 60% cut
Platelet – <b>PLT</b> (cells/l, cells/ $\mu$ l)	Number of thrombocytes (platelets) $PLT = PLT_{cal} \times (\text{cells/l, cells}/\mu\text{l})$
Mean Platelet Volume – <b>MPV</b> (fl)	Average volume of individual platelets derived from the PLT histogram
Thrombocrit – <b>PCT</b> (percentage, absolute)	Calculated from the PLT and MPV values $PCT_{percentage} = PLT \times MPV \times 100$ $PCT_{absolute} = PLT \times MPV$
White blood cell differential: <b>LYM, LY%</b> : lymphocytes <b>MID, MID%</b> : monocytes and some eosinophils <b>GRA, GR%</b> : neutrophil, eosinophil and basophil granulocytes	Absolute values counted in the channels determined by the three WBC discriminators:  Percentages calculated from the absolute WBC value.
Eosinophil – <b>EOS, EOS%</b> <b>EOS</b> shown by the second peak on the histogram	

## 4.4 Absolute and Linearity Ranges of Measured Parameters

The instrument is guaranteed to provide specified accuracy within its linearity range.

Beyond this linearity range, the instrument is able to display results, but may not guarantee accuracy characteristics.

If the value is over the maximum range of guaranteed linearity, the instrument cannot measure it and the result will be marked with an E (Error) flag.

To measure a sample, whose parameters exceed the maximum value indicated in the table below, pre-dilution is recommended. See section 5.3.1.3 of this manual.

The linearity ranges of primary parameters in normal measuring mode:

Parameter	Linearity Ranges	Maximum	Unit
WBC	0...100	150	10 <sup>9</sup> cells/liter
RBC	0...15	20	10 <sup>12</sup> cells/liter
PLT	0...700	1000	10 <sup>9</sup> cells/liter
HGB	0...250	400	g/l
HCT	0...100	-	%
MCV	30...150	-	fl
MPV	3...30	-	fl

**Table 2. Linearity ranges of parameters**

The linearity ranges for 1:5 pre-dilution mode:

Parameter	Linearity Ranges	Maximum	Unit
WBC	2...200	300	10 <sup>9</sup> cells/liter
RBC	1...30	40	10 <sup>12</sup> cells/liter
PLT	100...2000	3000	10 <sup>9</sup> cells/liter

**Table 3. Linearity ranges of Pre-dilution mode**

## 5. ROUTINE UTILIZATION and MEASURE

### 5.1 Sample handling

Since some time will usually elapse between collection of samples and counting, it is necessary to preserve the sample with an anti-coagulant to prevent large groups of cells forming into clots or lumps of cell matter that will clog the cell counter. Choice of anti-coagulant is very important, as some anticoagulants will affect the shape and size of blood cells. In general **EDTA**, preferably potassium based, is the only anti-coagulant recommended for use with electronic blood counters.

Care must be taken when using homemade containers pre-dosed with EDTA. If the container is not filled with enough blood, the ratio of EDTA to blood may reach a level, which results in osmotic transfer from the RBCs which shrinks them. **The ratio of EDTA to blood should not exceed 3 mg/ml.** Generally, we suggest using pre-manufactured sample tubes containing the necessary amount of EDTA. Also, when taking blood, please make sure that requirements attached to sample tubes are met.

---

**Important! Sample tubes must be filled to at least 7-8 mm height with blood otherwise correct sampling is not guaranteed!**

---

There is another possibility that can help the user in correct sampling by using the needle setting function. This is available in Measurement Local menu, and controls sampling height of the needle inside the sample tube. So, if you have a sample tube with a higher/lower bottom, you can control the sampling height adjusting this option. Needle offset is displayed in the lower left corner of the measurement screen.


<p><b><i>Attention! If you hurt yourself during analysis, biohazard substances can cause infection! Always use rubber gloves!</i></b></p>
---

To initiate analysis:

1. Invert the closed sample tube 11 times to achieve a homogenous sample.


*Do not shake the sample, because micro-bubbles can form inside which may cause erroneous sampling!*

*You have the possibility to use 3 different interchangeable adapters for different tube types. Tube types are shown in the next figures.*




**Figure 8. Tubes used in big adapter**

*Below you can see 3 types of tubes (microtainers) used in micro adapters. These are only examples given by us, you can try to use other type of microtainers as well.*



**Figure 9. Tubes used in micro adapter**



**Figure 10. Tube used in small adapter**

**1. Remove the cap!! It is very important because the tip will not pierce the cap!**

2. Position the sample tube in the sample rotor.
3. Push START key.


The sample rotor will turn into inside of the instrument and the needle draws sample from the tube. The aspirating needle is retracted, while its outer surface is automatically rinsed with diluent. This insures a low carry-over between samples. After a few seconds the rotor turns out. Now you can remove the sample tube from the adapter of the sample rotor.

## 5.2 Sample analysis

### 5.2.1 Sample preparation

Use K-EDTA anti-coagulated fresh whole blood as sample. **Prior** to sampling, mix the sample gently by **inverting it 11 times**. **Do not shake** as this could damage the blood cells.





### 5.2.2 Modifying lyse quantity

The default lyse quantity can be adjusted by pressing  on the MEASURE screen. Another option is to modify the lyse quantity by  $\pm 0.1$  ml or  $\pm 0.2$  ml is available during analysis.

Press  $\uparrow$  to **increase** the lyse quantity (**+0.1/0.2 ml**) if the separation between lysed RBCs and WBC populations is poorly differentiated, resulting in increased WBC and LYM counts.

Press  $\downarrow$  to **decrease** the lyse amount (**-0.1/0.2 ml**) if the WBC histogram seems to be shrunk to the left, i.e. the different WBC populations are overlapped. This can inhibit proper separation of WBC populations.

Patient limits			
Dog			
Lyse	0.65 ml	WBC method	A
ReLyse	0.00 ml	RBC method	A
Lyse V	3.50 ml		
Discr	80 fL		
WBC	6.00 - 17.0	$10^9/L$	

The two important parameters influencing lysing are lysing time and lyse quantity. You cannot change the lysing time, as it is adjusted to the lysing reagent supplied by Diatron.

Above you could read that the lyse quantity can be adjusted either at patient limits, or right before measurement.




## 5.2.3 Sample information

The software allows the user to enter information for each sample that has been, or will be measured. If an external PC keyboard (PS2) is used, it must be connected to the instrument **before** turning the instrument on.

Two options exist for sample information entry:

- immediately before analysis
- in the Database menu

To enter sample information prior to sample analysis, press the Measurement/Analysis key, and press **[#9999]** button on the screen. The following screen appears:



Sample information

Sample ID: 2      Date: 14.06.2007 11:23



Patient ID:



Species: Dog

Doctor:

Name:

Birth: 00.00.0000 Sex: -

Use the keypad OK key and the  and  keys to move through the settings.

Use the  and  keys to change the settings, and the keypad or external keyboard to enter text or numbers.

Press the  to accept data, cancel with  button.

Begin analysis by pressing the Start button.


---

**WARNING! Do not reach inside the instrument, as the needle can injure you!**






## 5.2.4 Results

When analysis is complete, the following screen is displayed, including all measured and calculated parameters as well as the WBC, RBC, PLT and VET5 histograms. Results and histograms will be stored automatically in the memory.



Horse: 1		11.06.2007 16:28	
WBC	8.73	RBC	8.15
LYM	1.66	LY%	19.0
MON	0.67	MO%	7.7
NEU	5.99	NE%	68.6
EOS	0.34	EO%	3.9
BAS	0.07	BA%	0.8
HGB	139	RDWc	20.1
HCT	36.53	MCH	17.0
MCV	45	MCHC	380
PLT	157	PDWc	31.0
MPV	7.1	PCT	0.11
Flags: None			



Use  and  keys on the screen to move through the results and histograms.

The following table summarizes **warning flags** and give an explanation of their possible cause and a few hints to overcome the problem:

**Uppercase letters refer to WBC or HGB problems:**

Flag	Meaning	Recommended user action
<b>E</b>	No WBC 3-part differential	Possible lyse problem. May occur in pathological lymphocytosis.
<b>H</b>	HGB blank is high, or no HGB blank	Repeat the blank measurement. If HGB blank is not stable there are probably bubbles in the WBC chamber: Run a cleaning and try blank again. Close the side door if open during measurement.
<b>B</b>	WBC blank is high, or no WBC blank	Repeat the blank measurement, or run prime lyse and try blank again. Possible lyse contamination, or noise problem.
<b>C, Q</b>	WBC clogging	Aperture clogging. Perform cleaning and repeat the measurement. If it is a general problem, please contact your Service Personnel. Low temperature reagents can cause it as well (mainly diluent), in this case you will have to wait until they reach room temperature.

**Table 4. Summary of warning flags related to WBC/HGB**

**Warning flags in lowercase refer to RBC or PLT problems:**

Flag	Meaning	Recommended user action
<b>p</b>	PLT blank is high, or no PLT blank	Run cleaning and repeat the blank measurement. Diluent or system cleanliness problem. If it is stable high, replace the diluent by opening a new tank.
<b>b</b>	RBC blank is high, or no RBC blank	Same action as in case of warning flag <b>p</b> .
<b>c</b>	RBC/PLT clogging	The same action as in case of the <b>C</b> warning flag.

**Table 5. Summary of warning flags related to EOS problems:**


Result Warning Flags		
Flag	Meaning	Description
<b>x</b>	<b>EOS Blank high</b>	Similar to the WBC B flag. Run a cleaning cycle and rerun blank.
<b>y</b>	<b>Clogging</b>	The Lyse2 measurement has encountered a clogging in the aperture. Run a cleaning cycle and repeat sample measurement
<b>z</b>	<b>Noise in the EOS Channel</b>	Possible lysing problem. Check if there is enough reagent in the container. Prime the Lyse2 reagent.

Warning flags can be grouped according to measurement conditions and according to the problems relating to the blood sample.

Measurement conditions: when the flags are related to clogging (**c**, **y**, **C**), probably hemolysing problems (**E**, **b**, **B**, **p**, **z**) and pressure problems (Fatal pressure error). In this case we suggest repeating the measurement.

The asterisk flag (\*) near a parameter shows some doubt suspected during the evaluation of that parameter. The reasons can be: a high PLT blank (PLT value will be marked), a case of indefinite discriminator setting (default location is used for some reasons, related parameters will be marked), etc.

Another flagging method is evaluation against the normal ranges. If some of the parameters is out of range it gets a (-) flag if under the range, or gets (+) if over the range. (And the given parameter will be highlighted as well.) You can customize ranges for all kind of patients by setting the corresponding lower and upper ranges. If you set 0 for a range limit, it will be not verified.

Horse: 5		11.06.2007 14:19	
WBC	7.47	RBC	6.48-
LYM	3.92	LY%	52.5
MON	0.04	MO%	0.5
NEU	3.40	NE%	45.5
EOS	0.10-	EO%	1.3
BAS	0.01	BA%	0.2
HGB	126	RDWc	20.1
HCT	30.37-	MCH	19.5
MCV	47	MCHC	416+
PLT	107	PDWc	33.1
MPV	7.3	PCT	0.08
Flags: None			
			

## 5.3 Measure


### 5.3.1 Measure local menu




### 5.3.2 Blank Measurement


Blank measurement is used for checking the cleanliness of the system and the reagents.

Blank measurement must be performed:

- Once daily, before sample analysis (this is done automatically before the first analysis in MEASURE function).
- After any reagent change (activated manually from the MEASURE/ MEASURE BLANK menu).
- After the replacement of any hardware component that is closely related to the measuring process (aspiration, dilution, counting, rinsing).

Press the Measure key and press the  button on the screen and select Measure blank.



When the Blank measurement is OK, press  to accept the result. The **Abacus junior vet 5** is then ready for analysis.

The instrument then displays a sample measurement screen, as shown, and is now ready to perform an analysis.

Cat: (5)		14.06.2007 16:31	
WBC		RBC	
LYM	LY%	HGB	RDWc
MON	MO%	HCT	MCH
NEU	NE%	MCV	MCHC
EOS	EO%	PLT	PDWc
BAS	BA%	MPV	PCT
Flags: None			
<div> <div></div> <div></div> <div></div> <div>#9999</div> <div></div> <div></div> </div>			

Cat: (5)		14.06.2007 16:32	
WBC		VETS/Probe:	
RBC		PLT	
		Lyse:	
Flags: None			
<div><div></div><div></div><div></div><div>#9999</div><div></div><div></div></div>			

There are 3 regions for blank value handling:

1. *Optimal* - all results are within acceptable ranges.
2. *Blank is high* - \* flag is displayed at relevant results.
3. *Blank exceeds acceptability* - no results displayed.

Parameter	1. No flag at parameter	2. * flag at result	3. E (error) flag at result
HGB	0-10 g/l	10 - 25 g/l	> 25 g/l
WBC	0 - 0.5 x10 <sup>3</sup> cells/μl	0.5 - 1.0 x10 <sup>3</sup> cells/μl	> 1.0 x10 <sup>3</sup> cells/μl
PLT	0 - 25 x10 <sup>3</sup> cells/μl	25 - 50 x10 <sup>3</sup> cells/μl	> 50 x10 <sup>3</sup> cells/μl
RBC	0 - 0.05 x10 <sup>6</sup> cells/μl	0.05 - 0.5 x10 <sup>6</sup> cells/μl	> 0.5 x10 <sup>6</sup> cells/μl
EOS	0-0.1 x10 <sup>3</sup> cells/μl	0.1-0.2 x10 <sup>3</sup> cells/μl	>0.2 x10 <sup>3</sup> cells/μl

**Table 6. Blank measurement ranges**

Accepted blank values are essential for proper calibration and quality control measurement. For this reason, no calibration or QC measurement can be performed without accepted blank values.

Quality control measurement and calibration can be performed **only** if all blank values are in the first region (receiving no flags or errors).

If analysis errors occur or the blank measurement is too high, an E error flag appears along with the affected parameter and “---” is displayed instead of results. In this situation perform a cleaning.



UNSUCCESSFUL BLANK MEASURE			
WBC	0.00	RBC	---E
LYM		HGB	0
MON		HCT	
NEU		MCV	
EOS	0.00	PLT	---E
BAS		MPV	
Flags: cbpCY			
<div> <div></div> <div></div> <div></div> <div></div> </div>			



UNSUCCESSFUL BLANK MEASURE			
WBC		VETS	
		Probe: R:938-93 W:938-93 D:937-93	
RBC		PLT	
		Lyse: W:0.50 D:2.00	
Flags: cbpCY			

### 5.3.3 Using Pre-diluted Mode


Perform an external pre-dilution of the sample using **Abacus junior vet 5** reagent diluent, or an isotonic saline solution. Dilute the sample to 1:5 ratio (1 part sample to 5 part diluent). Mix it well.

To perform the analysis:



1. Press the Measurement/Analysis key  then press the  soft key.



2. Select **Prediluted Mode**



3. Press  soft key. Notice that **Predil mode** now appears in lower left corner.



4. If you are processing a new sample. Press the **#9999** soft key and enter the patient information, then press **✓**.


Sample information	
Sample ID: 5	Date: 14.06.2007 16:34
Patient ID:	
Species:	Cat
Doctor:	
Name:	
Birth: 00.00.0000	Sex: -
<input type="button" value="X"/> <input type="button" value="✓"/>	

5. If you are processing a sample that has already been tested as non-prediluted: select **Repeat Last Sample** to avoid having to re-enter patient information.

Measure local menu	
1.	<b>Repeat last sample</b>
2.	Measure blank
<input type="checkbox"/> 3.	Prediluted mode
4.	Needle Height Setting - B: 0 mm
<input type="button" value="✓"/>	


6. Press Start button to begin the analysis. **Abacus junior vet 5** will automatically calculate the results with the 1:5 pre-dilution factor.


Horse: 1		11.06.2007 16:28	
WBC	8.73	Predil Mode	
LYM	1.66	LY%	19.0
MON	0.67	MO%	7.7
NEU	5.99	NE%	68.6
EOS	0.34	EO%	3.9
BAS	0.07	BA%	0.8
RBC	8.15	HGB	139
HCT	36.53	RDWc	20.1
MCV	45	MCH	17.0
PLT	157	MCHC	380
MPV	7.1	PDWc	31.0
PCT	0.11		
Flags: None			
<input type="button" value="←"/> <input type="button" value="⏏"/> <input type="button" value="#9999"/> <input type="button" value="→"/>			




## 6. DATABASE

Patient results are stored in the memory in chronological order, and can be retrieved at any time. Memorizing capacity is 2,000 measurements, including the complete parameter list, histograms, flags, sample data, and date/time of measurements. If memory is full, latest (actual) record will overwrite oldest record.

To access the Database table, press the Database key  on the front panel. The first screen that appears shows the most recent saved tests.



PatID	SmpID	Date
□8313	19	12.06.2007 13:40
□5612	20	12.06.2007 14:06
□3013	21	12.06.2007 14:16
□5630	22	12.06.2007 14:21
□5627	23	12.06.2007 14:32
□5457	24	12.06.2007 14:38

 #9999


### DATABASE

Pressing the ← or → key accesses the remaining, non-visible parameter results. The ↑ or ↓ key scrolls between the results individually.

To display the WBC, RBC, PLT and VET5 histograms, press the  soft key.

To print an individual result, highlight the result and press Print .

Database local menu	
1. Go to specified record	
2. Selection	▶
3. Change sort order	▶
4. Manage selected records	▶
5. View external	
6. Backup one day	



### <DATABASE LOCAL MENU>

From the database table screen, enter the Database local menu by pressing the **Menu** key.

The menu contains the following items:




Command	Function
<b>Go to specified record</b>	Jumps to a particular sample record. Enter the date and time, sample ID, and patient ID of the sample you want to view, and press <input checked="" type="checkbox"/> . The first sample meeting your parameters is then displayed. If you leave any ID blank (0), records are searched by date/time only
<b>Selection</b>	<b>Selects all</b> sample records in memory, or all having a specific date, time and ID. <b>Select by date, time and ID</b> allows you to select a range of specific records, and <b>Deselects all</b> deselects all records. Entering 0 as an ID searches by date/time only. Corresponding results are marked with a filled box.
<b>Change sort order</b>	Changes the order in which results are displayed: by time, sample ID or patient ID.
<b>Manage selected records</b>	<b>Sends selected</b> records to a PC, deletes them, or saves them to a diskette or USB device. Before <b>selecting Backup selected</b> records, insert a 3.5 inch diskette or connect a USB storage device. An empty floppy diskette can store data for 800 samples.
<b>View external</b>	Views previously saved data from a diskette or USB storage device.
<b>Backup one day</b>	Backs up all records from a specified day to a diskette or USB storage device. Select a day to backup, then press <input checked="" type="checkbox"/> to confirm.

## 7. UTILITIES


### 7.1 Maintenance

By selecting item (1) of the UTILITIES you can access the MAINTENANCE menu.



#### 7.1.1 Regular Maintenance Jobs

From Maintenance submenu, the user can initiate maintenance procedures such as cleaning, priming, draining chamber, reagent status or bleaching.




##### MAINTENANCE (1)

Select the required submenu.

##### 7.1.1.1 Cleaning


Item 1 in the above menu brings up cleaning functions.



Item 1 starts a washing cycle using the system cleaner reagent. This action is recommended if clogging problems are experienced (C or Q error flag).

Item 2 initiates a process that uses a light solution of hypochlorite (NaOCL), and washes the entire system with it. The instrument will ask for the cleaning solution in a sampling tube.

### 7.1.1.2 Priming



During the priming cycle, the fluidic system is rinsed with a large amount of diluent. It differs from the process in a start-up procedure; as in the latter case a simple filling up of the fluidics is performed. If fluid sensors are on, then the analyzer makes these procedures automatically, otherwise the User must initiate them activating the appropriate item within this submenu.



### 7.1.1.3 Draining chambers

Draining of chambers should be run before removal or replacement of parts related to the measuring chambers or apertures.

### 7.1.1.4 Reagent status

Reagent status					
Reagent pack installed: 01.01.2004/Reagent pack					
Samples: 0			Blanks: 0		
Diluent	Lyse	Cleaner	L_VETS	Rinse	Waste
9000	300	300	300	500	11000
100%	100%	100%	100%	100%	0%

The screen on the left shows reagent volumes in containers, as calculated by the instrument. As measurements are performed, the volumes are changing accordingly. When reagent volume in a container is running low, instrument will notify user, and ask replacement.

Press the  soft key, then press  to confirm.

The **Abacus junior vet 5** updates the installation date, reagent lifetime, and the amount of reagent in each container.

**WASTE HANDLING – VERY IMPORTANT**

*Waste contains poisonous substances (because of possible cyanide content) and human origin substances representing biohazard. These substances are representing potential danger to environment. For this reason, safe handling of the waste liquid is very important.*

*Please contact your distributor which kind of reagent is supplied to you, whether the lyse reagent contains cyanide or it is cyanide-free. Please, disregard point 1. below if your lyse reagent does not contain cyanide.*

*The following steps should be made for environmental protection and safety reasons:*

**1. Neutralization of waste containing cyanide:**

- *Use the set of cyanide neutralizer reagents (contact your distributor):*  
     *Component A (alkaline oxidant reagent)*  
     *Component B (neutralizer reagent)*
- *Put 5 ml/l of component A into the waste container when it is empty, then connect it to the analyzer.*
- *When the waste container becomes full replace it.*
- *Wait 1-2 hours for oxidation to complete.*
- *Put 5 ml/l of component B into the waste container.*
- *Wait 2-3 hours.*

**2. Neutralization of biohazard effect**


- *Independently of the fact that the waste contains cyanide or not, you should make this step.*
- *Put 2 ml/l hypochlorite solution into the waste. Close the cap, shake the container and wait 1 hour.*
- *Dispose of waste by spilling it into the drain system.*

**7.1.2 Weekly maintenance**


Weekly maintenance should be performed before turning on the power switch.

- **How to open the side doors:**

On the left side and the rear side of both instruments there is a side door, which gives access to the fluidic system and the mechanical parts easily (Figure 11 and 12). Other parts of the analyzer (electronic parts, etc.) can be accessed by opening the front cover and the rear cover.




**Figure 11-12.**



### ***Cleaning the washing head***

The washing head cleans the outer surface of the aspirating tip with saline diluent. Any salt build-up on the lower surface may cause malfunction during operation. Use a soft cloth dampened with water to clean this area. You can see the washing head indicated in the following figure:



***Figure13. Parts of measuring block***

1. Exit the Measurement menu. Open the side door after the needle has stopped moving.
2. Gently rub the lower surface of the washing head to remove the salt build-up.
3. Close the side door.

### ***Peristaltic pump maintenance***

The pump installed in the instrument is maintenance free. However, would you experience leakage from the pump, or vacuum error, you may replace the tube used.


Remove the tubes from the pump by opening the screw connectors.

Push in the two sides of the pump cassette.

Pull the cassette off the pump.

Pushing aside the two tubes, bend the tube out of the pump.

Now you have the tube as figure 1 shows.




Using a standard screwdriver, ease the tube on the plastic connector, and pull the tube off. Repeat the same procedure on the other end of the tube.

Now you can pull the plastic holders off the tube. Retain all parts except the old tube.

Using the new pump tube, slide the two grey plastic parts onto the new tube. Make sure they are aligned as shown on figure 1. Insert the 2 white plastic connectors into the ends of the tube. The new set must look like the tube on figure 1.

Put the tube back into the pump mechanics, and drive the grey plastic parts into their seats as figure 2 indicates. (view from the “top” of the pump)



Slide the cassette housing back. There are pins on the cassette, make sure to drive them in their paths.

The cassette must click into its place.

## 7.2 Calibration

The **Abacus junior vet 5** is pre-programmed to monitor DiatroCont 5 control. Performing QC determinations regularly verifies continued optimal performance.

It is recommended to do calibration in the following cases:

1. At analyzer installation, before beginning the analyses.
2. After replacing any component, related to the process of dilution or measurement.
3. When quality control measurements show any systematic error (bias) or they are outside predefined limits.
4. At regular time intervals (determined by the lab itself).
5. If you want to use the instrument in Prediluted mode.

Calibration can be performed in two ways:

1. User can enter calibration factors - without any calibration measurements – using the numerical keypad.
2. One-, two- or three-fold measurements of control or special calibrations material with known parameters. In this case, the instrument automatically calculates new factors using the following formula:

$$\text{New factor} = \frac{\text{Assigned value} \times \text{Stored factor}}{\text{Measured value(s) (or average of those)}}$$


---

**CAUTION!** New calibration will invalidate the previous factors. Old values cannot be retrieved, but can be reviewed in the VIEW CALIBRATIONS menu.

---

Calibration can be initiated by choosing Calibration in the UTILITIES.

1. Select the calibrator:






2. Enter the assay values for each parameter from the calibrator package insert.

Calibration

RBC	1.00
MCV	1.00
RDWc	1.00
PLT	1.00
MPV	1.00

Calibration

HGB	1.00
WBC	1.00
EOS	1.00
BAS	1.00

3. When all parameters are set, press the  soft key.

The display shows **Calib:(1/1)** at top left.

Calib:(1/1) 18.06.2007 10:53

WBC	LY%	RBC	
LYM	MO%	HGB	RDWc
MON	NO%	HCT	MCH
NEU	NE%	MCV	MCHC
EOS	EO%	PLT	PDWc
BAS	BA%	MPV	PCT

Flags: None

Calib:(1/1) 18.06.2007 10:53


WBC	VETSProbe:
RBC	PLT
	Lase:



Flags: None

Target values for calibrated parameters can be set within the following ranges:



Parameter	Low limit	High limit
RBC	1.00	8.00
HCT	0.1	0.6
MCV	50	120
RDW CV	10	50
PLT	30	800
PCT	0	2
MPV	5	15
PDW CV	5	50
HGB g/l	30	300
WBC	1.0	30.0

**Table 7. Calibration ranges**

9. When analysis is complete and the display shows the results, press  .



**Abacus junior vet 5** calculates and displays the new calibration factors.




## 7.3 Quality Control Procedure

By analyzing control materials, day-to-day reproducibility can be monitored. In this submenu, both target values and acceptable ranges for each parameter can be specified for different QC levels.

**NOTE:** Target values of the control material should be set only once, at the beginning of the QC measurements. Resetting parameters deletes previous QC results of the active level.

**CAUTION!** Any change in the QC material setting deletes previous QC results. It is strongly recommended to print results prior to changes.

1. In the Quality Control menu select **Set QC reference ranges**:



2. QC Reference values:

QC reference values

QC Type: 1

LOT No.: 0703-012N

Exp. date: 02.08.2007

Lyse 0.70 ml

WBC 7.30 +- 0.80 10<sup>9</sup>/l

RBC 4.46 +- 0.20 10<sup>12</sup>/l

3. Enter the value ranges for each test parameter as indicated on the paperwork accompanying the QC lot number you are using.

- use the keypad to change displayed values.
- to disable QC of a parameter, set it to 0.0
- use the PGDN and PGUP to view additional parameters.


4. Press  to accept the data, then press  to confirm.

5. Select Run QC to start the analysis. Press the Start button.

QC: QC 1/11 08.06.2007 11:25

WBC 7.17	LY% 3.06	LY% 3.06	RBC 4.49-	HGB 126+	RDWc 16.3+
LYM 3.06	MON 0.49	MO% 3.62	HCT 38.68-	MCH 28.1	MCHC 327
NEU 3.62	EOS 3.77+	EO% 0.11	MCV 86	PDWc 31.2	PCT 0.35
BAS 0.11	BA% 0.11		PLT 258		
			MPV 13.4		

Flags: None



Press  to accept the results and save them in the QC database

### Viewing the Accepted QC Database



1. Press the Utilities key, then select Quality Control.

2. Select **View table of QC measures.**

PatID	SmplID	Date	
QC 1/06	QC 1/06	08.06.2007	11:06
QC 1/07	QC 1/07	08.06.2007	11:10
QC 1/08	QC 1/08	08.06.2007	11:13
QC 1/09	QC 1/09	08.06.2007	11:17
QC 1/10	QC 1/10	08.06.2007	11:21
QC 1/11	QC 1/11	08.06.2007	11:25

## 7.4 Diagnostics

### Devices




### Self Test

The Self test is a procedure to verify proper operation of essential components of the instrument.

The Self test should be performed:


- At installation.
- After replacing any component.
- After extended time out of use.

1. Press the **Utilities** key and select **Diagnostics**.
2. Select **Self test**. **Abacus junior vet 5** then lists and checks subsystems.
3. When the test is finished, the **Abacus junior vet 5** displays a summary of the results.



## 7.5 Settings

Selecting item four (4) of the UTILITIES you can access this menu.



## 7.5.1 Printer settings


The following table lists the available printer selections, along with the printer language and specific printer models corresponding to each.

**Table 3-1: Compatible Printers**

Printer Selection	Printer Language	Supported Printer Models
Selko LPTH245 Built In	Special printer language	Built-in thermal printer
Epson ESC/P Raster	Epson ESC/P Raster	Epson Stylus series (select models only)
Canon BJC in 9-pin mode Canon BJC	Canon BJC	Canon BJ and BJC series (select models only)
HP DeskJet	PCL4	HP DeskJet series (select models only)
HP LaserJet	PCL4	HP LaserJet series and compatibles (LJ1100, etc.)
Epson 24-pin in 9-pin mode, Epson 24-pin	ESC/P2	Epson 24-pin dot matrix series (LQ-580, Epson 440, etc.)
Epson 9-pin	ESC/P	Epson 9-pin dot matrix series (FX-980, etc.)

To set up the instrument for your printer, go to the “**Utilities/Printer/Printer Settings**” menu. Select from the options using the up and down arrow keys within the text fields, and fill in the numerical fields using the number keys.

The general characteristics of the printable area of printer paper are below:





distance between reports.

The paper is defined by its size: it can be a standard size (A4, Letter) or any custom sized paper (specify the actual size). Printers cannot print on the whole surface of the paper. The blank area is described by the physical **Margins**, which may vary by printer models. The paper area inside the physical margins is called the *printable area*.


**Top margin** and **Left margin** settings are used for determining the location of printed results on the page. If more than one result is to be printed per sheet, use the **Vertical spacing** to specify the

On the first page of the Printer Settings sub-menu, printer type can be selected. Here you can choose between Seiko built-in printer, or any compatible external printer. Initially, only the printer driver can be selected, and when accepted, driver and printout format details become available.



Select the **Printer** matching your printer hardware. Pressing the  key will bring up printout details dialogs.

- ☐ Printout format: **Narrow full** (with histograms), or **Narrow text only** (without histograms)
- ☐ Table format: **Normal** or **Narrow**
- ☐ Printer cable: **COM2**
- ☐ Paper: **RollPaper**
- ☐ Units: **inch** or **cm**
- ☐ Paper size: **2.24 x 50.00**
- ☐ Left margin: **0.00**
- ☐ Top margin: **0.00**
- ☐ Printing mode: **Normal**



- Physical margin: **normal**
- One result per page: **No**
- Rollpaper: **Yes**
- Vertical spacing: **1.27**
- Autoprint: **Yes**
- Limits format: **L Graph H**

Printer settings	
Physical margin:	Normal
One result per page:	No
Rollpaper:	Yes
Vertical spacing:	1.27
Autoprint:	No
Limits format:	L Graph H

- Print flags: **Yes**
- Print warning flags: **Yes**
- Clogging report: **No**
- Serial number on result: **Yes**

Printer settings	
Print flags:	Yes
Print warning flags:	Yes
Clogging report:	No
Serial number on result:	Yes

- For an external printer, the following default settings are recommended:
  - ☐ Printout format: **Full with histograms**, or **Wide text only** (results are printed without histograms)
  - ☐ Paper: select as appropriate for your printer (**Note:** Never select **Rollpaper** for an HP LaserJet.)
  - ☐ Left and top margins: **1.00**
  - ☐ Rollpaper: **No** (especially for HP LaserJets)

Configure any other settings as needed. See Table 3-2 for a complete list.

**Table 3-2: Printer Settings**

Setting	Description
<b>Printer</b>	Printer used with the VET5 — see Table 3-1.
<b>Printout format</b>	Overall format used for printouts. Available selections depend on the selected <b>Printer</b> . (See page 3-6 for examples.)
<b>Table format</b>	Printed width of tables included in the printouts.
<b>Printer cable</b>	Interface used for an external printer: LPT (parallel) or USB.
<b>Paper</b>	Paper type used in the selected <b>Printer</b> .
<b>Units</b> <b>Paper size</b> <b>Left margin</b> <b>Top margin</b>	Default units of measure, paper size, and print margins for the selected <b>Paper</b> . These can be customized if you're using a non-standard paper. Left and top margins control the position of printed results on the page.

**Table 3-2: Printer Settings**

Setting	Description
<b>Printing mode</b>	Overall printout size ( <b>Mini</b> , <b>Small</b> , <b>Normal</b> , <b>Enhanced</b> , or <b>Large</b> ), width ( <b>Wide</b> or <b>Narrow</b> ), and print speed ( <b>Fast</b> ). The available modes depend on the selected <b>Printer</b> .
<b>Physical margin</b>	Margin correction for printers that have smaller printable areas than usual. Use <b>Normal</b> unless the right edge of the printed text is missing or appears in the next line.
<b>One result per page</b>	Select <b>Yes</b> to print each analysis on a separate page.
<b>Rollpaper</b>	Select <b>Yes</b> if using the built-in printer.
<b>Vertical spacing</b>	The vertical distance between printed results when multiple reports are printed on each page.
<b>Autoprint</b>	Select <b>Yes</b> to print all results automatically after each analysis.
<b>Limits format</b>	The format for printing the reference ranges: <ul style="list-style-type: none"> <li>To print numeric values, select the displayed range (for example, [4.51...5.64]).</li> <li><b>None</b> prints no range values.</li> <li><b>Graphical</b> prints a graph showing where the results falls in range.</li> <li><b>L-H Graph</b> prints the low and high limits of the range, followed by a graph showing where the result falls in the range.</li> <li><b>Graph L-H</b> prints a graph showing where the result falls in the range, followed by the low and high limits of the range.</li> <li><b>L Graph H</b> prints the low limit of the range, a graph showing where the result falls in the range, and the high limit of the range.</li> </ul>
<b>Print flags</b>	Select <b>Yes</b> to include in the printout any measured values that fell outside the reference range, as well as any errors during analysis.
<b>Print warning flags</b>	Select <b>Yes</b> to include in the printout all flags that appeared in the results.
<b>Clogging report</b>	Select <b>Yes</b> to include information on aperture clogging on all printouts.
<b>Serial number on result</b>	Select <b>Yes</b> to include the VET5's serial number on all printouts.

Press the  soft key to accept the settings.




## 7.5.2 Customize


Use the Customize menu to set the language and date format used by the **Abacus junior vet 5**, along with the screen saver delay, and the date range used for printing combined results from the **Abacus junior vet 5**.

1. Press the **Utilities** key. 


2. Select **Settings**



3. Select **Customize**



4. Select **General Settings**



This is a collection of settings influencing instrument operation, customization.




Table 3-4 describes the available settings.

**Table 3-4: General Settings**

Screen saver wait time	0 to 30 minutes (0 = screen saver off)
Language	English. International customers: please contact your distributor for information about available languages.
PC Link	Offline, USB, Baud rate (150 to 115200). Protocol used by the <b>VET5</b> to communicate with a connected computer. – If the <b>VET5</b> is not connected to a computer, choose Offline. – If the <b>VET5</b> is connected to a computer using a USB cable, choose USB.
Serial Protocol Version	1.0, 1.7, 2.20, 2.23, 3.0, or 3.1. Protocol versions used by practice management software to interface with the .
Serial's Decimal Dot	Dot (period), slash, comma, or lang def ( as defined by language)


VSx Link	Offline, USB, RS232. Protocol used by the <b>Vet5</b> to communicate with a connected VetScan Classic or VS2 for combining results. – If the <b>Vet5</b> is not connected to a VetScan Classic or VS2, choose Offline. – If the <b>Vet5</b> is connected to a VetScan VS2 through a USB cable, choose USB. – If the <b>Vet5</b> is connected to VetScan Classic through a serial cable, choose RS232
Date format of VetScan	When printing combined <b>Vet5</b> and VetScan results, the date format used by the VetScan: MM/DD/YY or DD/MM/YY (DD = date, MM = month, YY = year). Set both instruments to the same format.

## 7.5.3 Units


Use the unit menu to set the measurement units the **Abacus junior vet 5** will use.

1. Press the Utilities key 


2. Select **Settings**.



3. Select **Customize**.



4. Select **Units**.



5. Set the units as needed:

- Use the ◀ and ▶ keys to select a unit.
- Use the ▲ and ▼ keys to select a setting for that unit.

Table 3-5 lists the available unit settings.

**Table 3-5: Units**


Count unit	cells/liter ( <b>cells/l</b> ), cells/μl ( <b>cells/μl</b> ). "Prec." displays a more precise value.
HGB unit	grams/liter ( <b>g/l</b> ), grams/deciliter ( <b>g/dl</b> ), millimols/liter ( <b>mmol/l</b> )
PCT, HCT unit	percentage ( <b>%</b> ), absolute ( <b>abs</b> )
RDW, PDW mode	standard deviation ( <b>sd</b> ), coefficient of variation ( <b>cv</b> )

6. Press the ☒ on-screen key to save the settings.

## Laboratory information settings

Enter your clinic or laboratory information as follows, to be printed automatically on report headers. You can also use this procedure to edit or change this information.

Enter the name and address of your clinic or laboratory. Press OK to move to the next line. You can enter up to 40 characters in each line.




Press ☒ to accept the settings.


## User mode settings

If you want to require users to log in and enter passwords to use **Abacus junior vet 5**, and to be able to track individual usage of the unit, you can enable the **Abacus junior vet 5**'s Multi user mode. This will require each user to have a unique user ID and password.

1. Press the **Utilities** key.
2. Select **Settings**.
3. Select **Customize**.
4. Select **User mode**.




5. Follow the on-screen instructions to enable multi user mode, to add or modify users, or to change a user's password.





## 7.5.4 Date and time


The date and time of each analysis is stored with the results. This menu allows setting the built-in clock and the format of the date displayed.


1. Press the **Utilities** key, and then select **Settings**.
2. Select **Date and time**.
3. Select Set date and time.




4. Type in the date and time, then press .



5. Select formats for displaying the date (item 2, 3 or 4) and time (item 5 or 6), then press .




6. Press  to accept the settings.

## 8. Printing

This chapter covers information on making reports on measured samples.

### 8.1 Printouts

When required, the following items can be sent to an external printer or to a built-in printer by pressing the  function key button.

- \* Database result(s) (table format)
  - \* Database (specified patient results with histograms)
  - \* QC result (Levey-Jennings chart)
  - \* QC result(s) (table format)
  - \* Calibration results
  - \* Last measured blank result
  - \* Last measured patient result (with histograms)
  - \* Last measured QC result
  - \* Device information and statistics
  - \* Self test result
  - \* Set parameters

The appropriate printout format can be selected in UTILITIES/SETTINGS/PRINTER SETTINGS).

## Thermal paper printout

ABACUS JUNIOR VET 5

Patient ID:  
3590

Name:

Mode:  
Horse

Doctor:

/ -

Birth/Sex: 00.00.0000

Report date: 20.06.2007 S/N: 260005

Test time: 12.06.2007 16:45

Test Result Sample ID: 37

WBC	9.03	$10^9/l$	5.4	14.3
LYM	2.68	$10^9/l$	1.5	7.7
MON	0.47	$10^9/l$		1.5
NEU	5.61	$10^9/l$	2.3	9.5
EOS	0.25	$10^9/l$	0.1	1
BAS	0.03	$10^9/l$		0.5
LY%	29.7	%	17	68
MO%	5.2	%		14
NE%	62.1	%	22	80
EO%	2.8	%	1	8
BA%	0.3	%		3
RBC	9.88	$10^{12}/l$	6.8	12.9
HGB	171	g/l	110	190
HCT	42.54	%	32	53
MCV	43	fL	37	59
MCH	17.3	pg	12.3	19.7
MCHC	401	g/l	310	390
RDWc	22.3	%		
PLT	79	$10^9/l$	100	400
PCT	0.05	%		
MPV	6.3	fL		
PDWc	32.6	%		

WBC


EOS

RBC

PLT


## Letter printout

### ABACUS JUNIOR VET 5









## 9. Fluidic Schematics



## 10. Warning labels on the analyzer

Label	Meaning	Explanation	Chapter
	<b>Biohazard</b>	The sample and the waste are potentially infectious material	2.2.4
	<b>Corrosive</b>	Reagents may cause corrosion or skin irritation.	2.2.4; 2.3
	<b>Warning</b>	General warning of injury.	
	<b>Sharp needle warning</b>	The sampling needle may cause injury.	5.2.3

**Table 6. Warning labels**